

PMManager Release Notes

Version 4.00 (.08) Released 12-May-2026

New features:

- x940 communication support

Updates:

- Support for X940
- Updated "Newport User Commands"
- Updated build files to use Visual Studio 2022
- MSVC v14 Redistributable that comes with VS 2022
- Minor internal updates

Fixes:

- Minor internal Fixes

Known issues:

- X940 Ethernet not yet operational
- Some HW settings in x940 / x938 are not fully controllable by PMManager
 - Trigger settings
 - TTL, Analog out (Raw/Digital), Analog Filter
 - Range Type (Physical / Logical)
- Time-out when searching for non-existing devices on Bluetooth or Ethernet can take up to a minute!
- COM-Object bug mainly affecting Python applications

Version 3.93 (.02) Released 23-December-2024

New features:

- Supports upgrading 819-WL sensor firmware
Note that 819-WL sensors are not otherwise supported by this version

Updates:

- Add auto-ranging option for line graphs

Fixes:

-

Known issues:

-

Version 3.92 (.03) Released 11-July-2024

New features:

- Added 818-RAD support for X938 on PMManager
- Enhancement: Add auto scale option for line graphs
- Extended "Time Frame" of graphs to be up to 99:59 hours

Updates:

- Update "Newport User Commands"
- COM objects version 10.7

Fixes:

- PMManager previous version didn't run on Win 7
- When attempting to open a second instance of PMManager there was an assertion failure

Known issues:

-