




Model	Dimensions [in. (mm)]		Thread
	A	B	C
AJS100-05K	1.41 in. (35.8 mm)	0.50 in. (12.7 mm)	1/4-100
AJS127-0.5K	1.41 in. (35.8 mm)	0.50 in. (12.7 mm)	1/4-127
AJS100-.75K	1.85 in. (47.0 mm)	0.75 in. (19.1 mm)	1/4-100
AJS100-1K	2.27 in. (57.5 mm)	1.00 in. (25.4 mm)	1/4-100