





1/4-20 x 1" S[ACING, 1/2" BORDER FOR VIS-
M6 x 25mm SPACING, 12.5mm BORDER FOR M-VIS



NOTE:
CASTER HEIGHT ADJUSTMENT: ± 0.2 "
ADDITIONAL HEIGHT ADJUSTMENT FOR ISOLATOR VERTICAL STROKE: ± 0.3 "

VIS SYSTEM	DIM A	DIM B	DIM C in.	DIM D in.
VIS3030-XXX-I	30.0 in	30.0 in	32.8	32.8
VIS3036-XXX-I	30.0 in	36.0 in	32.8	38.8
VIS3048-XXX-I	30.0 in	48.0 in	32.8	50.8
VIS3636-XXX-I	36.0 in	36.0 in	38.8	38.8
VIS3648-XXX-I	36.0 in	48.0 in	38.8	50.8
VIS3660-XXX-I	36.0 in	60.0 in	38.8	62.8
VIS3672-XXX-I	36.0 in	72.0 in	38.8	74.8
M-VIS3030-XXX-I	750 mm	750 mm	32.8	32.8
M-VIS3036-XXX-I	750 mm	900 mm	32.8	38.8
M-VIS3048-XXX-I	750 mm	1200 mm	32.8	50.8
M-VIS3636-XXX-I	900 mm	900 mm	38.8	38.8
M-VIS3648-XXX-I	900 mm	1200 mm	38.8	50.8
M-VIS3660-XXX-I	900 mm	1500 mm	38.8	62.8
M-VIS3672-XXX-I	900 mm	1800 mm	38.8	74.8

XXX = SG2, SG4, PG2, OR PG4