

StarViewer for Android Release Notes

Version 3.10.2 Released on 26-March-2026

New features:

- Restored Juno Support via USB (previously broken with Android updates)
- Ariel & StarViewer log files formatted to be compatible with StarLab

Updates:

- Adjusted screen partitioning after android-15 update changes
- Updated with latest Android APK libraries (31, 33 & 35)
- Other minor updates & fixes

Fixes:

- Truncated lines in log
- App doesn't close properly
- Quasar battery display
- Watts & Joules measurements confusion
- Fix for measurements streaming (broken after some time)

Known issues:

-