

1.1.2.13 Accessories for Fan Cooled Thermal Sensors

1.1.2.13.1 Protective Covers with Target Pattern for F150A-BB-26, FL250A-BB-50 & FL400A-BB-50 Sensors

All the protective covers are made of black anodized aluminum, and have a cross pattern for alignment. These protective covers are for F150A-BB-26, FL250A-BB-50 & FL400A-BB-50 sensors and are not supplied with the sensors. They can be ordered separately for these sensors.

The covers protect the sensor absorbers when the sensors are not in use, and during the preparations for using the sensor they allow easy alignment with the laser beam.

Sensor F150A-BB-26
with Protective Cover




Sensor FL250A-BB-50 / FL400A-BB-50
with Protective Cover



Protective Cover	For use with	Part Number
F150A-26 Protective Cover	F150A-26 sensor	1G05448
FL250A-50 & FL400A-50 Protective Cover	FL250A-BB-50-V1, FL250A-BB-50-V1 with 5m cable, FL400A-BB-50, FL400A-BB-50 with 5m cable sensors	1G05449

F150A-26 Protective Cover



FL250A-50 & FL400A-50 Protective Cover

