

## StarLab Release Notes

---

### Version 4.00 (.08) Released 12-May-2026

#### *New features:*

- Support for automatic conversion of Ariel logfiles uploaded by StarLab or StarViewer, and native logfiles from StarViewer. Reads and displays the following log file types:
  - StarLab native log
  - StarLab-uploaded Ariel log
  - StarViewer/Android native log
  - StarViewer/Android-uploaded Ariel log
  - StarViewer/iOS-uploaded Ariel log

#### *Updates:*

- Updated "Ophir User Commands"
- Updated build files to use Visual Studio 2022
- MSVC v14 Redistributable that comes with VS 2022
- Minor internal updates

#### *Fixes:*

- Log upload progress bar fixed
- Minor internal fixes

#### *Known Issues:*

- Time-out when searching for non-existing devices on Bluetooth or Ethernet can take up to a minute!

## Version 3.93 (.02) Released 23-December-2024

### *New features:*

- The COM object now supports wireless connections to Ariel

### *Updates:*

- COM objects version 10.8
- Updated "Ophir User Commands"

### *Fixes:*

- Allow DOSAGE mode for StarLite, with PD300-RM

### *Known Issues:*

- Time-out when searching for non-existing device on Bluetooth or Ethernet can take up to a minute!

## Version 3.92 (.03) Released 11-July-2024

### *New features:*

- Added PD300RM support for Centauri on StarLab
- Enhancement: Add auto scale option for line graphs
- Support connecting to Centauri via Ethernet
- Extended "Time Frame" of graphs to be up to 99:59 hours

### *Updates:*

- Update "Ophir User Commands"
- COM objects version 10.7

### *Fixes:*

- StarLab previous version didn't run on Win 7
- When attempting to open a second instance of StarLab there was an assertion failure

### *Known Issues:*

- Time-out when searching for non-existing device on Bluetooth or Ethernet can take up to a minute!